


# Electrode Code Chart

## Subdural Electrodes

REF: FG40J-SP10X-040

Platinum  
Marker



**Contacts are shown not exposed (recording side DOWN).**




USA

Rx ONLY

|        |
|--------|
| RED    |
| ORANGE |
| YELLOW |
| GREEN  |
| BLUE   |
| VIOLET |
| WHITE  |
| BLACK  |
| GRAY   |
| BROWN  |

The subdural electrodes should be handled with care to prevent damage.  
(A direct pull or stress on the electrode may cause a loss of contact recordings).

**For Single Use Only. Do Not Re-Sterilize or Reuse. For Surgical Use Only.**  
**Not Intended for Implantation (≥ 30 days). Do not use if packaging is damaged.**



**AD-TECH** MEDICAL INSTRUMENT CORPORATION

1901 WILLIAM STREET, RACINE, WISCONSIN 53404

Phone: 262.634.1555

Email: sales@adtechmedical.com