


PMT[®] Corporation

Electrode Product Guide

Featuring:

- Fully color coded tails
- Platinum discs and wires
- Enhanced surface options on grid and strip electrodes for greater pliability
- Custom capabilities
- Local product specialists
- Variety of connectors for intraoperative and LTM



Committed to advancements in intracranial electrodes since 1979

PMT[®]
PMT CORPORATION

Quality, Innovation and Customer Driven

PMT[®] Corporation
1500 Park Road~Chanhassen~MN~55317~USA
PHONE: 1-952-470-0866~FAX:1-952.470.0865
Toll Free in USA Only: 1-800-MANKIND (626-5463)
www.pmtcorp.com, www.permark.com
Copyright © 2009 All Rights Reserved

Table of Contents PMT® Intracranial Electrodes:

We invite and thank you for previewing PMT® electrodes. All of our electrodes are crafted using the highest quality materials in a class 10,000 cleanroom. Please consult your local PMT® Product Specialist or distributor for inquiries of electrode variances not shown in the catalog.

<u>SUBDURAL AND DEPTH ELECTRODES</u>	<u>PAGE</u>
Subdural Strip and Grid Product Overview.....	3
Cortac® Subdural Strip Electrodes.....	4
Cortac® Subdural Grid Electrodes.....	5-9
Cortac® Subdural Strip and Grid Dual Sided Electrodes.....	10-11
Depthalon® Depth Electrodes and Anchor Bolts.....	12
Macro/Micro and Micro Array Electrodes.....	13

SPHENOID ELECTRODES AND CABLES

(KITS INCLUDE 2 EACH ELECTRODES AND NEEDLES WITH CRIMPED OR NON-CRIMPED OPTIONS)

Stainless Steel Sphenoidal® Electrodes With Cables.....	14
Stainless Steel Sphenoidal® Electrodes Without Cables.....	14
Platinum Sphenoidal® Electrodes With Cables.....	14
Platinum Sphenoidal® Electrodes Without Cables.....	14
Silver Sphenoidal® Electrodes With Cables.....	14
Silver Sphenoidal® Electrodes Without Cables.....	14
Sphenoidal® Electrode Cable 1.5 mm Socket Set.....	14
Crimping Tool for Non-Crimped Sphenoidal® Electrode.....	14

DEPTH PLACEMENT KITS AND ACCESSORIES

Leksell Stereotactic System® Frames.....	15
(SPLIT CANNULA, T-DRIVER, SPLIT CANNULA BUSHINGS, T-DRIVER BUSHINGS)	
Medtronic StealthStations® Stereotactic Frames.....	15
(SPLIT CANNULA, T-DRIVER)	
BRW™/CRW™ Stereotactic Frames.....	16
(SPLIT CANNULA, T-DRIVER, SPLIT CANNULA BUSHING, DRILL GUIDE BUSHINGS, T-DRIVER BUSHINGS)	
Electrode Trocar and Passing Needles.....	16

CUSTOM ORDERING FORMS

Subdural Custom Ordering Form.....	18
Depth Custom Ordering Form.....	19
Cortac Microwire Ordering Form.....	20

In 1984 PMT® pioneered the development and introduction of reliable, precision Cortac® Cortical Strip and Grid Electrodes, Depthalon® Depth Electrodes and Sphenoidal® Electrodes used in EEG monitoring. These products are used to evoke potential stimulation and electrocorticography diagnostic applications that are critical for determining appropriate surgical approaches to treating epilepsy patients. PMT® Corporation's commitment to innovative products and customer satisfaction has made PMT® a leader in epilepsy monitoring, intraoperative and LTM. PMT platinum cortical surface electrodes, and microwire electrodes are approved by the FDA for LTM up to 30 days. Utilizing the finest quality materials in the manufacturing process, PMT® GUARANTEES to deliver you the quality and service you need.

PMT® Corporation is proud to be a CE Marked, FDA and EN 46001/ISO 13485 registered manufacturer. Federal Law (U.S.A.) restricts the devices shown in this catalog for sale by or on the order of a physician. All PMT products are covered by one of the best comprehensive product liability insurances available to the medical manufacturing industry. Contact your local PMT® Product Specialist or distributor for more information.

All electrodes are supplied sterile for single use only. Drawings in the catalog are not to scale and dimensions listed are approximate. Model numbers are subject to change without notice. U.S. Patent 6162101, 6024702, 5961909, 5902236 and other patents pending.


Committed to advancements in intracranial electrodes since 1979

PMT®
PMT CORPORATION


Quality, Innovation and Customer Driven

PMT® Corporation
1500 Park Road~Chanhassen~MN~55317~USA
PHONE: 1-952-470-0866~FAX:1-952.470.0865
Toll Free in USA Only: 1-800-MANKIND (626-5463)
www.pmtcorp.com, www.permark.com
Copyright © 2009 All Rights Reserved

PMT[®] Strip and Grid Electrode Product Overview



Parts of the Strip or Grid



I) The grids are available in single or multiple exits so the electrode may be “trimmed” if desired.

II) Electrode discs are available in 2.0 mm, 3.0 mm, or 4.5 mm size.

III) The electrode disc may also be darkened or radiopaque for x-rays.


IV) This electrode has 8 columns and 8 rows and can also be described as an 8x8.

V) Electrode discs and wiring are available in 100% (Series 2110) platinum or stainless steel (Series 2111) options.

VI) Strips are available in .5 mm, .76 mm or 1.01 mm thickness. Grids are available in .5 mm or .76 mm thickness.

VII) Grids and strips can be ordered with original or enhanced surface (TX). Enhanced surface offers greater pliability which is particularly useful on grids.

Parts of the Tail




I) The “tail contacts” portion of the tail connects to the PMT[®] In-Line Cable (Series 2125).

PMT[®] Cortac[®] Strip Electrodes:

Unless otherwise specified, default electrode disc size is 3.0 mm with a .76 mm strip thickness, original surface, and random color coded tails. For other options, please consult your PMT[®] Representative or Distributor.


Model Number	Description*	Material	Cable Model Number
2110-02-032	2 contact, 1x2, one exit	Platinum	2125-02-01 In-Line
2111-02-032	2 contact, 1x2, one exit	Stainless Steel	2125-02-01 In-Line




Model Number	Description* †	Material	Cable Model Number
2110-04-032	4 contact, 1x4, one exit	Platinum	2125-04-01 In-Line
2111-04-032	4 contact, 1x4, one exit	Stainless Steel	2125-04-01 In-Line



Model Number	Description*	Material	Cable Model Number
2110-06-032	6 contact, 1x6, one exit	Platinum	2125-06-01 In-Line
2111-06-032	6 contact, 1x6, one exit	Stainless Steel	2125-06-01 In-Line




Model Number	Description*	Material	Cable Model Number
2110-08-032	8 contact, 1x8, one exit	Platinum	2125-08-01 In-Line
2111-08-032	8 contact, 1x8, one exit	Stainless Steel	2125-08-01 In-Line




PMT[®] Cortac[®] Grid Electrodes: Unless otherwise specified, default electrode disc size is 3.0 mm with a .5 mm grid thickness, enhanced surface and random color coded tails. For other options, please consult your PMT[®] Representative or Distributor.

Model Number	Description*	Material	Cable Model Number
2110TX-10-074	10 contact, 2x5, one exit	Platinum	2125-10-01 In-Line
2111TX-10-074	10 contact, 2x5, one exit	Stainless Steel	2125-10-01 In-Line




Model Number	Description*	Material	Cable Model Number
2110TX-12-076	12 contact, 2x6, one exit	Platinum	2125-12-01 In-Line
2111TX-12-076	12 contact, 2x6, one exit	Stainless Steel	2125-12-01 In-Line



Model Number	Description†	Material	Cable Model Number
2110TX-12-090	12 contact, 2x6, two exits	Platinum	2125-06-02 In-Line
2111TX-12-090	12 contact, 2x6, two exits	Stainless Steel	2125-06-02 In-Line



Model Number	Description*	Material	Cable Model Number
2110TX-16-076	16 contact, 2x8, one exit	Platinum	2125-16-01 In-Line
2111TX-16-076	16 contact, 2x8, one exit	Stainless Steel	2125-16-01 In-Line




PMT[®]
PMT CORPORATION

Quality, Innovation and Customer Driven


PMT[®] Corporation
1500 Park Road~Chanhassen~MN~55317~USA
PHONE: 1-952-470-0866~FAX:1-952.470.0865
Toll Free in USA Only: 1-800-MANKIND (626-5463)
www.pmtcorp.com, www.permark.com
Copyright © 2009 All Rights Reserved

PMT[®] Cortac[®] Grid Electrodes: Unless otherwise specified, default electrode disc size is 3.0 mm with a .5 mm grid thickness, enhanced surface and random color coded tails. For other options, please consult your PMT[®] Representative or Distributor.


Model Number	Description	Material	Cable Model Number
2110TX-16-074	16 contact, 4x4, one exit	Platinum	2125-16-01 In-Line
2111TX-16-074	16 contact, 4x4, one exit	Stainless Steel	2125-16-01 In-Line



Model Number	Description	Material	Cable Model Number
2110TX-20-074	20 contact, 4x5, one exit	Platinum	2125-10-02 In-Line
2111TX-20-074	20 contact, 4x5, one exit	Stainless Steel	2125-10-02 In-Line



Model Number	Description*	Material	Cable Model Number
2110TX-32-074	32 contact, 4x8, one exit	Platinum	2125-16-02 In-Line
2111TX-32-074	32 contact, 4x8, one exit	Stainless Steel	2125-16-02 In-Line



PMT[®]
PMT CORPORATION


Quality, Innovation and Customer Driven

PMT[®] Corporation
1500 Park Road~Chanhassen~MN~55317~USA
PHONE: 1-952-470-0866~FAX:1-952.470.0865
Toll Free in USA Only: 1-800-MANKIND (626-5463)
www.pmtcorp.com, www.permark.com
Copyright © 2009 All Rights Reserved


PMT[®] Cortac[®] Grid Electrodes:

Unless otherwise specified, default electrode disc size is 3.0 mm with a .5 mm grid thickness, enhanced surface and random color coded tails. For other options, please consult your PMT[®] Representative or Distributor.


Model Number	Description	Material	Cable Model Number
2110TX-32-079	32 contact, 4x8, two exits	Platinum	2125-16-02 In-Line
2111TX-32-079	32 contact, 4x8, two exits	Stainless Steel	2125-16-02 In-Line



Model Number	Description†	Material	Cable Model Number
2110TX-32-090	32 contact, 4x8, four exits	Platinum	2125-08-04 In-Line
2111TX-32-090	32 contact, 4x8, four exits	Stainless Steel	2125-08-04 In-Line




Model Number	Description* ‡	Material	Cable Model Number
2110TX-64-051	64 contact, 8x8, one exit	Platinum	2125-16-04 In-Line
2111TX-64-051	64 contact, 8x8, one exit	Stainless Steel	2125-16-04 In-Line




PMT[®] Cortac[®] Grid Electrodes:

Unless otherwise specified, default electrode disc size is 3.0 mm with a .5 mm grid thickness, enhanced surface and random color coded tails. For other options, please consult your PMT[®] Representative or Distributor.

Model Number	Description*	Material	Cable Model Number
2110TX-64-074	64 contact, 8x8, two exits	Platinum	2125-16-04 In-Line
2111TX-64-074	64 contact, 8x8, two exits	Stainless Steel	2125-16-04 In-Line



Model Number	Description	Material	Cable Model Number
2110TX-64-079	64 contact, 8x8, four exits	Platinum	2125-16-04 In-Line
2111TX-64-079	64 contact, 8x8, four exits	Stainless Steel	2125-16-04 In-Line




PMT[®] Cortac[®] Grid Electrodes:

Unless otherwise specified, default electrode disc size is 3.0 mm with a .5 mm grid thickness, enhanced surface and random color coded tails. For other options, please consult your PMT[®] Representative or Distributor.

Model Number	Description ‡	Material	Cable Model Number
2110TX-64-090	64 contact, 8x8, eight exits	Platinum	2125-08-08 In-Line
2111TX-64-090	64 contact, 8x8, eight exits	Stainless Steel	2125-08-08 In-Line


‡ Grid Thickness is .76 mm




PMT[®] Dual Sided Strip and Grid Electrodes:

Unless otherwise specified all dual strips and grid thickness are 1 mm. For other options, please consult your PMT[®] Representative or Distributor.


Model Number	Description	Material	Cable Model Number
2110-12-020	12 contact, dual sided grid, 4.5 mm electrode disc two exits <i>mini connection</i>	Platinum	2120-06-02 Mini Cable
2111-12-020	12 contact, dual sided grid, 4.5 mm electrode disc two exits <i>mini connection</i>	Stainless Steel	2120-06-02 Mini Cable



Model Number	Description	Material	Cable Model Number
2110-16-020	16 contact, dual sided grid, 4.5 mm electrode disc two exits <i>mini connection</i>	Platinum	2120-08-02 Mini Cable
2111-16-020	16 contact, dual sided grid, 4.5 mm electrode disc two exits <i>mini connection</i>	Stainless Steel	2120-08-02 Mini Cable



Model Number	Description	Material	Cable Model Number
2110-24-019	16 contact, dual sided grid, 3.0 mm electrode disc two exits <i>mini connection</i>	Platinum	2120-12-02 Mini Cable
2111-24-019	16 contact, dual sided grid, 3.0 mm electrode disc two exits <i>mini connection</i>	Stainless Steel	2120-12-02 Mini Cable
2110-24-020	16 contact, dual sided grid, 4.5 mm electrode disc two exits <i>mini connection</i>	Platinum	2120-12-02 Mini Cable
2111-24-020	16 contact, dual sided grid, 4.5 mm electrode disc two exits <i>mini connection</i>	Stainless Steel	2120-12-02 Mini Cable



PMT[®]
PMT CORPORATION


Quality, Innovation and Customer Driven

PMT[®] Corporation
1500 Park Road~Chanhassen~MN~55317~USA
PHONE: 1-952-470-0866~FAX:1-952.470.0865
Toll Free in USA Only: 1-800-MANKIND (626-5463)
www.pmtcorp.com, www.permark.com
Copyright © 2009 All Rights Reserved

PMT[®] Dual Sided Strip and Grid Electrodes:

Unless otherwise specified all dual strips and grid thickness are 1 mm. For other options, please consult your PMT[®] Representative or Distributor.

Model Number	Description	Material	Cable Model Number
2110-48-088	48 contact, dual sided grid, 3.0 mm electrode disc one exit	Platinum	2125-12-04 In-Line
2111-48-088	48 contact, dual sided grid, 3.0 mm electrode disc one exit	Stainless Steel	2125-12-04 In-Line




PMT[®] Depthalon[®] Electrodes: Unless otherwise specified, a stainless steel stylet is provided, (tungsten available upon request) default electrode disc size is 5.0 mm with 4 and 6 contact electrodes and 2.0 mm disc size with 8 and 10 contact electrodes. Contact “1” is located directly on the distal end of the electrode. For other options, please consult your PMT[®] Representative or Distributor.

Model Number	Description	Material	Cable Model Number
2102-04-IL-SS	4 contact	Platinum	2125-04-01 In-Line
2101-04-IL-SS	4 contact	Stainless Steel	2125-04-01 In-Line




Model Number	Description	Material	Cable Model Number
2102-06-IL-SS	6 contact	Platinum	2125-06-01 In-Line
2101-06-IL-SS	6 contact	Stainless Steel	2125-06-01 In-Line



Model Number	Description	Material	Cable Model Number
2102-08-IL-SS	8 contact	Platinum	2125-08-01 In-Line
2101-08-IL-SS	8 contact	Stainless Steel	2125-08-01 In-Line



Model Number	Description	Material	Cable Model Number
2102-10-IL-SS	10 contact	Platinum	2125-10-01 In-Line
2101-10-IL-SS	10 contact	Stainless Steel	2125-10-01 In-Line



PMT[®] Anchor bolts: Made of titanium with four lengths and two burr hole options. Supplied sterile in double tyvek pouches.

Model Number	Name	Description of Use
2103-23	Anchor Bolts (3.2 mm Burr Hole)	Secures the depthalon electrode to the skull (25 mm, 30 mm, 35 mm, 40 mm)
2103-27	Anchor Bolts (2.7 mm Burr Hole)	Secures the depthalon electrode to the skull (25 mm, 30 mm, 35 mm, 40 mm)

Model Number 2103-23



Model Number 2103-27




PMT[®]
PMT CORPORATION

Quality, Innovation and Customer Driven

PMT[®] Corporation
1500 Park Road~Chanhassen~MN~55317~USA
PHONE: 1-952-470-0866~FAX:1-952.470.0865
Toll Free in USA Only: 1-800-MANKIND (626-5463)
www.pmtcorp.com, www.permark.com
Copyright © 2009 All Rights Reserved


PMT[®] Macro/Micro Array Electrodes: Unless otherwise specified, macro micro array electrodes have .76 mm grid thickness, 3.0 mm electrode disc diameter and platinum material. **Other configurations available. Please consult your PMT[®] Representative or Distributor for more information.**

Model Number	Description	Material	Cable Model Number
2112-24-128-2	24 Macro electrode discs, 128 micro electrodes, <i>size 50</i> micron	Platinum	Not Applicable
2112-24-128-3	24 macro electrode discs, 128 micro electrodes, <i>size 75</i> micron	Platinum	Not Applicable



PMT[®] Micro Array Electrodes: Unless otherwise specified, micro array electrode contacts are built with platinum material. **Other configurations available. Please consult your PMT[®] Representative or Distributor for more information.**

Model Number	Description	Material	Cable Model Number
2112-00-16-2	16 micro electrodes, <i>size 50</i> micron, in-line	Platinum	2125-16-1 In-Line
2112-00-16-3	16 micro electrodes, <i>size 75</i> micron, in-line	Platinum	2125-16-1 In-Line



PMT[®]
PMT CORPORATION

Quality, Innovation and Customer Driven


PMT[®] Corporation
1500 Park Road~Chanhassen~MN~55317~USA
PHONE: 1-952-470-0866~FAX:1-952.470.0865
Toll Free in USA Only: 1-800-MANKIND (626-5463)
www.pmtcorp.com, www.permark.com
Copyright © 2009 All Rights Reserved

PMT[®] Sphenoidal[®] Electrodes: Can be supplied individually or in sets, using 36 gauge platinum, stainless steel or silver wire. Unless specified otherwise sphenoidal electrodes arrive in 30 cm length, with 3 - 5 mm tip. For other options, please consult your PMT[®] Representative or distributor.

Model Number	Material Description	(With Cables)	Crimped
2101-31-02	2 Each Stainless Steel Electrodes, 2 Needles	WITH 2 Cables	Yes
2101-32-02	2 Each Stainless Steel Electrodes, 2 Needles	WITH 2 Cables	No
2102-31-02	2 Each Platinum Electrodes, 2 Needles	WITH 2 Cables	Yes
2102-32-02	2 Each Platinum Electrodes, 2 Needles	WITH 2 Cables	No
2103-31-02	2 Each Silver Electrodes, 2 Needles	WITH 2 Cables	Yes
2103-32-02	2 Each Silver Electrodes, 2 Needles	WITH 2 Cables	No

Model Number	Material Description	(Without Cables)	Crimped
2101-31-03	2 Each Stainless Steel Electrodes and 2 Needles		Yes
2101-32-03	2 Each Stainless Steel Electrodes and 2 Needles		No
2102-31-03	2 Each Platinum Electrodes and 2 Needles		Yes
2102-32-03	2 Each Platinum Electrodes and 2 Needles		No
2103-31-03	2 Each Silver Electrodes and 2 Needles		Yes
2103-32-03	2 Each Silver Electrodes and 2 Needles		No

Model Number	Material Description
2122	Sphenoidal Electrode EEG Cable 1.5mm socket Set (2 cables)
2105	Crimping Tool for Non-Crimped Sphenoidal Electrode



PMT[®] Depthalon[®] Depth Placement Kits: The depth placement kits are supplied non-sterile. For other options, please consult your PMT[®] Representative or distributor. Drawings are not to scale.

Depth Electrode Placement Kit for Leksell Stereotactic System[®] Frames. This Kit Includes:

<u>Model Number</u>	<u>Name</u>	<u>Description of Use</u>
2103-39	Split Cannula	Guides the depthalon electrode to the anchor bolt
2103-26	T-Driver	Used to screw in the anchor bolt
2103-62	Upper Split Cannula Bushing	Holds the split cannula in the stereotactic system
2103-61	Lower Split Cannula Bushing	Holds the split cannula in the stereotactic system
2103-58	Upper Guide Bushing for T-Driver	Guides the T-Driver in the stereotactic system
2103-59	Lower Guide Bushing for T-Driver	Guides the T-Driver in the stereotactic system

Model Number 2103-39




Model Number 2103-26



Model Number 2103-58




Model Number 2103-59



Model Number 2103-62



Model Number 2103-61



Depth Electrode Placement Kit for Medtronic StealthStations[®] Stereotactic Frames. This Kit Includes:

<u>Model Number</u>	<u>Name</u>	<u>Description of Use</u>
2103-40	Split Cannula	Guides the depthalon electrode to the anchor bolt
2103-38	T-Driver	Used to screw in the anchor bolt

Model Number 2103-40



Model Number 2103-38



PMT[®]
PMT CORPORATION

Quality, Innovation and Customer Driven

PMT[®] Corporation
1500 Park Road~Chanhassen~MN~55317~USA
PHONE: 1-952-470-0866~FAX:1-952.470.0865
Toll Free in USA Only: 1-800-MANKIND (626-5463)
www.pmtcorp.com, www.permark.com
Copyright © 2009 All Rights Reserved

PMT[®] Depthalon[®] Depth Placement Kits: The depth placement kits are supplied non-sterile. For other options, please consult your PMT[®] Representative or distributor. Drawings are not to scale.

Depth Electrode Placement Kit for BRW/CRW[™] Stereotactic Frames


This Kit Includes:

<u>Model Number</u>	<u>Name</u>	<u>Description of Use</u>
2103-40	Split Cannula	Guides the depthalon electrode to the anchor bolt
2103-26	T-Driver	Used to screw in the anchor bolt
2103-47	Split Cannula Bushing	Holds the split cannula in the stereotactic system
2103-48	Drill Guide Bushing	Guides the drill in the stereotactic system
2103-45	Guide Bushing for T-Driver	Guides the T-Driver in the stereotactic system
2103-51	Thumbscrew	Used to secure the split cannula


Model Number 2103-40



Model Number 2103-26




Model Number 2103-47




Model Number 2103-51



Model Number 2103-48



Model Number 2103-45




PMT[®] Passing Needle and Trocar: Used to tunnel electrode tails away from the surgical site

Trocar for Electrodes with Mini Connection Model 2103-55



In-line Passing Needle Model 2105-5



PMT[®]
PMT CORPORATION

Quality, Innovation and Customer Driven

PMT[®] Corporation
1500 Park Road~Chanhassen~MN~55317~USA
PHONE: 1-952-470-0866~FAX:1-952.470.0865
Toll Free in USA Only: 1-800-MANKIND (626-5463)
www.pmtcorp.com, www.permark.com
Copyright © 2009 All Rights Reserved

PMT[®] CORPORATION

Custom Ordering Forms

Please consult your local distributor or sales representative for additional styles not shown in the catalog or feel free to complete the custom ordering forms on the following pages. Fax or mail to us at...



Fax: +1 952.470.0865, Subject: Custom Order Request

Address:

PMT Corporation
Attn: Custom Order Request
1500 Park Road
Chanhassen MN 55317
United States

PMT[®]Cortac[®] Custom Ordering Form: Please complete the below chart according to your specifications. Drawings are not to scale.

Hospital: _____ Doctor: _____ Contact: _____
 Patient Name: _____ Sales Rep.: _____
 Doctor's Signature: _____ Date: _____ Order #: _____



Please draw contact configuration needed on grid.

1. Contacts

- Number of Contacts: _____
- Contact Diameter: 4.5mm 3mm 2mm
- Contact Material: Platinum Stainless Steel
- Center to Center Contact Spacing (DIM A): _____ mm ± _____ mm
- Contact Numbering: All Top & Bottom Rows (If other configuration is desired, please mark in figure above).

2. Dimensional

- Length of Grid / Strip (DIM B): _____ mm ± _____ mm
- Width of Grid / Strip (DIM C): _____ mm ± _____ mm
- Thickness of Grid / Strip: _____ mm ± _____ mm
- Length of Lead (DIM D): _____ mm ± _____ mm
- Number of Tails: _____
- Length of Tails (DIM E): _____ mm ± _____ mm

3. Connector:

- Connector Type: Standard Heavy Duty Mini In - Line
- Connector Markings Type: # ID Tail Color Coding

4. Radiopaque:

- Radiopaque Contacts: YES NO (If yes, mark contacts to be radiopaqued in figure above).
- Radiopaque Grid: YES NO


5. Miscellaneous:

Please identify any other specifications required: _____

PMT[®] Depthal[®] Depth Custom Ordering Form:

Please complete the below chart according to your specifications. Drawings are not to scale.

Hospital: _____ Doctor: _____ Contact: _____
 Patient Name: _____ Sales Rep.: _____
 Doctor's Signature: _____ Date: _____ Order #: _____



1. Number of Contacts: _____
2. Contact Material: _____ Platinum _____ Stainless Steel
3. Contact Length: _____ 5mm _____ 2mm
4. Center to Center Spacing of the Electrodes (DIM A): _____ CM
5. Length from Tip Electrode to Collar (DIM B): _____ CM
6. Length of Collar (DIM C): _____
7. Length from Collar to Connector (DIM D): _____
8. Main Tubing Length (DIM E): _____ CM
9. Connector Heat Shrink Color Coding: _____ YES _____ NO
10. Connector Type: _____ Standard _____ Mini _____ Heavy Duty _____ In-Line
11. Stylet Entry Position (DIM F): _____ CM
12. Stylet Material Preference: _____ Stainless Steel _____ Tungsten
13. Stylet Length (DIM G): _____ CM
14. Miscellaneous Specifications Needed: _____

FOR PMT CORPORATION USE ONLY

SHIPMENT VERIFICATION


Engineering: _____
 Shipping Manager: _____
 Quality Assurance Manager: _____

PMT[®] Cortac[®] Microwire Custom Ordering Form:

Please complete the below chart according to your specifications. Drawings are not to scale.

Hospital: _____ Doctor: _____ Contact: _____
 Patient Name: _____ Sales Rep.: _____
 Doctor's Signature: _____ Date: _____ Order #: _____

Please draw Micro/Macro
contact configuration needed.



1. Contacts

- Number of Microwire Contacts: _____
- Microwire Contact Configurations: _____ X
- Microwire Contact Diameter: 40um 75um 125um 175um
- Number of Macro Contacts: _____
- Macro Contact Diameter: 2mm 3mm 4.5mm
- Center to Center Microwire Contact Spacing (DIM A): _____ mm
- Center to Center Macro Contact Spacing (DIM A): _____ mm
- Contact Numbering: All Top & Bottom Rows (If other configuration is desired, please mark in figure above).

2. Dimensional

- Length of Grid / Strip (DIM B): _____ mm
- Width of Grid / Strip (DIM C): _____ mm
- Thickness of Grid / Strip: _____ mm
- Number of Leads: _____
- Length of Leads (DIM D): _____ mm
- Number of Tails: _____
- Length of Tails (DIM E): _____ mm

3. Radiopaque

- Radiopaque Contacts: YES NO
- Radiopaque Grid: YES NO

(If yes, mark contacts to be radiopaqued in figure above)

PMT[®]
PMT CORPORATION

Quality, Innovation and Customer Driven

PMT[®] Corporation
 1500 Park Road~Chanhassen~MN~55317~USA
 PHONE: 1-952-470-0866~FAX:1-952.470.0865
 Toll Free in USA Only: 1-800-MANKIND (626-5463)
 www.pmtcorp.com, www.permark.com
 Copyright © 2009 All Rights Reserved



PMT[®] WARRANTY AND DISCLAIMER:

PMT[®] Corporation warrants that reasonable care was used in the choice of materials and manufacture of this product. PMT[®] Corporation shall not be liable for any incidental or consequential loss, damage, expense, or injury directly or indirectly arising from the use of this product, other than the replacement of it, should PMT[®] Corporation's investigation show the product was defective at the time of shipment by PMT[®] Corporation. Because PMT[®] Corporation has no control over the conditions of use, patient selection, post-surgical stresses, or handling of the device after it leaves the possession of PMT[®] Corporation, PMT[®] Corporation cannot guarantee the outcome or effects of any procedure involving PMT[®] Neurosurgical Electrodes. The foregoing warranties, as conditioned and limited are in lieu of and exclude all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of MERCHANTABILITY OR FITNESS FOR USE. PMT[®] Corporation neither assumes nor authorizes any other person to assume for it any other or additional liability or responsibility in connection with this product.

Caution:

Federal Law (USA) restricts this device to sale by or on the order of a physician.

Committed to advancements in intracranial electrodes since 1979

PMT[®]
PMT CORPORATION

Quality, Innovation and Customer Driven

PMT[®] Corporation
1500 Park Road~Chanhassen~MN~55317~USA
PHONE: 1-952-470-0866~FAX:1-952.470.0865
Toll Free in USA Only: 1-800-MANKIND (626-5463)
www.pmtcorp.com, www.permark.com
Copyright © 2009 All Rights Reserved

PN 12876 Rev. 1 11/09